


Bottle[®]

CORONARY BIFURCATION ANGIOPLASTY CATHETER

with diameters corresponding
to anatomical dimensions
of the coronary bifurcation vessels

One balloon
with two diameters



BALTON
since
1980

BALTON[®]

CE₂₄₀₉

CORONARY BIFURCATION ANGIOPLASTY CATHETER

Angioplasty catheter Bottle is dedicated to:

- Bifurcation optimization after classic stent implantation.
- Angioplasty procedure in coronary bifurcation vessels.

Angioplasty catheter Bottle:

- One balloon with two diameters corresponding to anatomical dimensions of the coronary bifurcation vessel.
- Placement precision thanks to three radiopaque markers.

Technical data

Balloon material: nylon
Nominal pressure: 10 atm
Catheter coating: hydrophilic
Pushability and trackability: excellent
Inflation and deflation time: very short
Proximal shaft diameter: 2F (0.65 mm)
Distal shaft diameter: 2.7F (0.90 mm)
Compatible guidewire: max. .014"
Compatible guiding catheter: 5F (.056")
Catheter length: 140 cm
Delivery system type: Rapid Exchange

Caution:

This publication could include technical or other inaccuracies or typographical errors. Changes are periodically added to the information on herein; these changes will be incorporated in the new editions of the publication. Balton Ltd. may make improvements and/or changes in the product(s) described in this publication at any time. This publication is only a commercial information and cannot be handled the commercial offer in the meaning of the art.66 §1 of Polish Civil Code.

Reference:	Nominal balloon diameter [mm]		Nominal balloon length [mm]
	Dist	Prox	
CBB2.50x2.75x10140	2.50	2.75	10
CBB2.50x3.25x10140	2.50	3.25	10
CBB2.75x3.00x10140	2.75	3.00	10
CBB2.75x3.50x10140	2.75	3.50	10
CBB3.00x3.50x10140	3.00	3.50	10
CBB3.00x3.75x10140	3.00	3.75	10
CBB3.50x4.25x10140	3.50	4.25	10
CBB3.75x4.50x10140	3.75	4.50	10
CBB2.50x2.75x15140	2.50	2.75	15
CBB2.50x3.25x15140	2.50	3.25	15
CBB2.75x3.00x15140	2.75	3.00	15
CBB2.75x3.50x15140	2.75	3.50	15
CBB3.00x3.50x15140	3.00	3.50	15
CBB3.00x3.75x15140	3.00	3.75	15
CBB3.50x4.25x15140	3.50	4.25	15
CBB3.75x4.50x15140	3.75	4.50	15
CBB2.50x2.75x18140	2.50	2.75	18
CBB2.50x3.25x18140	2.50	3.25	18
CBB2.75x3.00x18140	2.75	3.00	18
CBB2.75x3.50x18140	2.75	3.50	18
CBB3.00x3.50x18140	3.00	3.50	18
CBB3.00x3.75x18140	3.00	3.75	18
CBB3.50x4.25x18140	3.50	4.25	18
CBB3.75x4.50x18140	3.75	4.50	18
CBB2.50x2.75x23140	2.50	2.75	23
CBB2.50x3.25x23140	2.50	3.25	23
CBB2.75x3.00x23140	2.75	3.00	23
CBB2.75x3.50x23140	2.75	3.50	23
CBB3.00x3.50x23140	3.00	3.50	23
CBB3.00x3.75x23140	3.00	3.75	23
CBB3.50x4.25x23140	3.50	4.25	23
CBB3.75x4.50x23140	3.75	4.50	23