

DIALYSIS

FULL RANGE OF CATHETERS FOR SHORT TERM DIALYSIS ACCESS
TO MEET EVERY ACCESS NEED



BALTON[®]

BALTON
since
1980


CE 0197

Dialysis set is intended for obtaining of access to a large blood vessel for dialysis performance.



Benefits of dialysis catheters:


- Fully radiopaque catheter introduced by the Seldinger technique
- Easy to place and remove
- Soft, radiopaque distal tip designed to reduce the likelihood of vessel trauma
- Rotatable security wings enhancing the patient comfort
- Specifications on catheter size and priming volume reported on the device
- Good elasticity and kinking resistant
- Sterile access design



Standard set components



Optional set components



DIALYSIS SET, SINGLE LUMEN, STRAIGHT


| Most commonly used sizes

Reference: guide wire "J" type	Reference: guide wire Straight type	Catheter size	Catheter length [cm]	Dilatator I/II	Guide wire size [inch]	Guide wire length [cm]
ZDJP8F10	ZDJP8F10P	8F	10	8F /10F	.035"	60
ZDJP8F15	ZDJP8F15P		15			
ZDJP8F20	ZDJP8F20P		20			

NOTE: On demand, it is possible to select additional components. While ordering, add the following endings to the reference index:
 S – scalpel; Y – Y hub; SR – Guiding syringe 5 ml; SL – Syringe 5 ml; KLLT – TEGO connector; KLLI – Injection cap; IY – V needle; N – Nitinol guide wire; H – Hydrophilic coating
 E.g.: ZDJP8F15HYKLLT.

| **On request** (please contact us for more information)

ZDJP8F (length from 10 to 30 cm)



DIALYSIS SET, SINGLE LUMEN "Y"

Reference: guide wire "J" type	Reference: guide wire Straight type	Catheter size	Catheter length [cm]	Dilatator I/II	Guide wire size [inch]	Guide wire length [cm]
ZDJY8F10	ZDJY8F10P	8F	10	8F /10F	.035"	60
ZDJY8F15	ZDJY8F15P		15			
ZDJY8F18	ZDJY8F18P		18			
ZDJY8F20	ZDJY8F20P		20			


NOTE: On demand, it is possible to select additional components. While ordering, add the following endings to the reference index:
 S – scalpel; Y – Y hub; SR – Guiding syringe 5 ml; SL – Syringe 5 ml; KLLT – TEGO connector; KLLI – Injection cap; IY – V needle; N – Nitinol guide wire; H – Hydrophilic coating
 E.g.: ZDJY8F20HPYKLLT.

DIALYSIS SET

Standard set components



Optional set components



DIALYSIS SET, DOUBLE LUMEN, STRAIGHT TYPE

Most commonly used sizes

Reference: guide wire "J" type	Reference: guide wire Straight type	Catheter size	Catheter length [cm]	Dilatator I/II	Guide wire size [inch]	Guide wire length [cm]	Version
ZDD6.5F10	ZDD6.5F10P	6.5F	10	7F / 8F	.027"	60	Pediatric
ZDD6.5F12	ZDD6.5F12P		12				
ZDD7F10	ZDD7F10P	7F	10				
ZDD7F15	ZDD7F15P		15				
ZDD8F10	ZDD8F10P	8F	10	8F / 10F	.035"	For Adults	
ZDD8F12	ZDD8F12P		12				
ZDD8F15	ZDD8F15P		15				
ZDD9F10	ZDD9F10P	9F	10				
ZDD9F15	ZDD9F15P		15				
ZDD11F15	ZDD11F15P	11F	15				10F / 12F
ZDD11F18	ZDD11F18P		18				
ZDD11F20	ZDD11F20P		20				
ZDD11.5F15	ZDD11.5F15P	11.5F	15				
ZDD11.5F20	ZDD11.5F20P		20				
ZDD12F15	ZDD12F15P	12F	15	12F/14F	.035"	For Adults	
ZDD12F18	ZDD12F18P		18				
ZDD12F20	ZDD12F20P		20				
ZDD13F15	ZDD13F15P	13F	15				
ZDD13F18	ZDD13F18P		18				
ZDD13F20	ZDD13F20P		20				

NOTE: On demand, it is possible to select additional components. While ordering, add the following endings to the reference index:
 S – scalpel; Y – Y hub; SR – Guiding syringe 5 ml; SL – Syringe 5 ml; KLLT – TEGO connector; KLLI – Injection cap; IY – V needle; N – Nitinol guide wire; H – Hydrophilic coating
 E.g.: ZDD12F20HPYKLLT.

On request (please contact us for more information)

ZDD6.5F (length from 7 to 20 cm)
 ZDD7F (length from 7 to 20 cm)
 ZDD8F (length from 7 to 20 cm)
 ZDD9F (length from 10 to 20 cm)

ZDD10F (length from 10 to 20 cm)
 ZDD11F (length from 12 to 45 cm)
 ZDD11.5F (length from 12 to 45 cm)
 ZDD12F (length from 12 to 45 cm)


ZDD13F (length from 12 to 45 cm)
 ZDD13.5F (length from 12 to 45 cm)
 ZDD14F (length from 12 to 45 cm)



Standard set components



Optional set components



DIALYSIS SET, DOUBLE LUMEN, CURVED TYPE

Most commonly used sizes

Reference: guide wire "J" type	Reference: guide wire Straight type	Catheter size	Catheter length [cm]	Dilatator I/II	Guide wire size [inch]	Guide wire length [cm]	Version
ZDD6.5F10Z	ZDD6.5F10ZP	6.5F	10	7F / 8F	.027"	60	Pediatric
ZDD6.5F12Z	ZDD6.5F12ZP		12				
ZDD8F10Z	ZDD8F10ZP	8F	10	8F / 10F	.035"		
ZDD8F12Z	ZDD8F12ZP		12				
ZDD8F15Z	ZDD8F15ZP		15				
ZDD11F15Z	ZDD11F15ZP	11F	15	10F / 12F	.035"	70	For Adults
ZDD11F18Z	ZDD11F18ZP		18				
ZDD11F20Z	ZDD11F20ZP		20				
ZDD12F15Z	ZDD12F15ZP	12F	15	10F / 12F			
ZDD12F18Z	ZDD12F18ZP		18				
ZDD12F20Z	ZDD12F20ZP		20				

NOTE: On demand, it is possible to select additional components. While ordering, add the following endings to the reference index:
S – scalpel; **Y** – Y hub; **SR** – Guiding syringe 5 ml; **SL** – Syringe 5 ml; **KLLT** – TEGO connector; **KLLI** – Injection cap; **IY** – V needle; **N** – Nitinol guide wire; **H** – Hydrophilic coating
 E.g.: ZDD12F20HZPSYKLLT.

On request (please contact us for more information)

ZDD6.5F_Z (length from 7 to 20 cm)
 ZDD7F_Z (length from 7 to 20 cm)
 ZDD8F_Z (length from 7 to 20 cm)
 ZDD9F_Z (length from 10 to 20 cm)

ZDD10F_Z (length from 10 to 20 cm)
 ZDD11F_Z (length from 12 to 45 cm)
 ZDD11.5F_Z (length from 12 to 45 cm)
 ZDD12F_Z (length from 12 to 45 cm)


ZDD13F_Z (length from 12 to 45 cm)
 ZDD13.5F_Z (length from 12 to 45 cm)
 ZDD14F_Z (length from 12 to 45 cm)

DIALYSIS SET

Standard set components



Optional set components



DIALYSIS SET, TRIPLE LUMEN, STRAIGHT TYPE


Most commonly used sizes

Reference:	Catheter size	Catheter length [cm]	Dilatator I/II	Guide wire size [inch]	Guide wire length [cm]
ZDT12F15	12F	15	10F /12F	.035"	70
ZDT12F18		18			
ZDT12F20		20			

NOTE: On demand, it is possible to select additional components. While ordering, add the following endings to the reference index:
 S – scalpel; Y – Y hub; SR – Guiding syringe 5 ml; SL – Syringe 5 ml; KLLT – TEGO connector; KLLI – Injection cap; IY – V needle; N – Nitinol guide wire; H – Hydrophilic coating
 E.g.: ZDT12F15HKLLIYN.

On request (please contact us for more information)

ZDT12F (length from 12 to 45)



DIALYSIS SET, TRIPLE LUMEN, CURVED TYPE

Most commonly used sizes

Reference:	Catheter size	Catheter length [cm]	Dilatator I/II	Guide wire size [inch]	Guide wire length [cm]
ZDT12F15Z	12F	15	10F /12F	.035"	70
ZDT12F18Z		18			
ZDT12F20Z		20			

NOTE: On demand, it is possible to select additional components. While ordering, add the following endings to the reference index:
 S – scalpel; Y – Y hub; SR – Guiding syringe 5 ml; SL – Syringe 5 ml; KLLT – TEGO connector; KLLI – Injection cap; IY – V needle; N – Nitinol guide wire; H – Hydrophilic coating
 E.g.: ZDT12F15HZPSYKLLT.

On request (please contact us for more information)

ZDT12F_Z (length from 12 to 45)

Caution:

The content of this document, including provided data, information, trademarks and logos are the sole property of BALTON. Any reproduction, whether in whole or in part, is strictly prohibited and will be considered a violation of BALTON's copyrights and other intellectual proprietary rights. ©2021 Balton Sp. z o.o. All rights reserved. This document and any images contained therein is intended solely for the use by healthcare professionals and BALTON's distributors.

This publication is for information purposes only and should not be construed as a commercial offer in the meaning of art. 66(1) of the Polish Civil Code.

All medical devices made commercially available by BALTON must be used in accordance with the enclosed instructions for use.

This publication presents only a small sample of our complete offering of medical devices. For further information please contact your Balton representative or email us at balton@balton.pl.



Bioline Comercial Ltda.
 R. José Getúlio, 579 14° andar cj 141/146
 Liberdade | SP | 01509 001
 +55 11 3274 7600
www.bioline.com.br