




Stryderyk

CORONARY ANGIOPLASTY CATHETER
WITH DELIVERY SYSTEM RAPID EXCHANGE (RX)



- Balloon material – polyamide
- Precision of manufacturing – easy to use
- High pressure – resistant
- Allows to dilate even the highly resistant lesion
- Optimalization of size, shape and durability
- Short inflation and deflation time




Reference:	Nominal balloon diameter [mm]	Nominal balloon length [mm]
CBA1.25x10140	1.25	10
CBA1.25x15140		15
CBA1.25x20140		20
CBA1.25x25140		25
CBA1.25x30140		30
CBA1.25x35140		35
CBA1.25x40140	40	
CBA1.5x10140	1.5	10
CBA1.5x15140		15
CBA1.5x20140		20
CBA1.5x25140		25
CBA1.5x30140		30
CBA1.5x35140		35
CBA1.5x40140	40	
CBA1.75x10140	1.75	10
CBA1.75x15140		15
CBA1.75x20140		20
CBA1.75x25140		25
CBA1.75x30140		30
CBA1.75x35140		35
CBA1.75x40140	40	
CBA2.0x10140	2.0	10
CBA2.0x15140		15
CBA2.0x20140		20
CBA2.0x25140		25
CBA2.0x30140		30
CBA2.0x35140		35
CBA2.0x40140	40	
CBA2.25x10140	2.25	10
CBA2.25x15140		15
CBA2.25x20140		20
CBA2.25x25140		25
CBA2.25x30140		30
CBA2.25x35140		35
CBA2.25x40140	40	
CBA2.5x10140	2.5	10
CBA2.5x15140		15
CBA2.5x20140		20
CBA2.5x25140		25
CBA2.5x30140		30
CBA2.5x35140		35
CBA2.5x40140	40	
CBA2.75x10140	2.75	10
CBA2.75x15140		15
CBA2.75x20140		20
CBA2.75x25140		25
CBA2.75x30140		30
CBA2.75x35140		35
CBA2.75x40140	40	
CBA3.0x10140	3.0	10
CBA3.0x15140		15
CBA3.0x20140		20
CBA3.0x25140		25
CBA3.0x30140		30
CBA3.0x35140		35
CBA3.0x40140	40	

Reference:	Nominal balloon diameter [mm]	Nominal balloon length [mm]
CBA3.25x10140	3.25	10
CBA3.25x15140		15
CBA3.25x20140		20
CBA3.25x25140		25
CBA3.25x30140		30
CBA3.25x35140		35
CBA3.25x40140	40	
CBA3.5x10140	3.5	10
CBA3.5x15140		15
CBA3.5x20140		20
CBA3.5x25140		25
CBA3.5x30140		30
CBA3.5x35140		35
CBA3.5x40140	40	
CBA3.75x10140	3.75	10
CBA3.75x15140		15
CBA3.75x20140		20
CBA3.75x25140		25
CBA3.75x30140		30
CBA3.75x35140		35
CBA3.75x40140	40	
CBA4.0x10140	4.0	10
CBA4.0x15140		15
CBA4.0x20140		20
CBA4.0x25140		25
CBA4.0x30140		30
CBA4.0x35140		35
CBA4.0x40140	40	
CBA4.5x10140	4.5	10
CBA4.5x15140		15
CBA4.5x20140		20
CBA4.5x25140		25
CBA4.5x30140		30
CBA4.5x35140		35
CBA4.5x40140	40	
CBA5.0x10140	5.0	10
CBA5.0x15140		15
CBA5.0x20140		20
CBA5.0x25140		25
CBA5.0x30140		30
CBA5.0x35140		35
CBA5.0x40140	40	
CBA6.0x10140	6.0	10
CBA6.0x15140		15
CBA6.0x20140		20
CBA6.0x25140		25
CBA6.0x30140		30
CBA6.0x35140		35
CBA6.0x40140	40	
CBA7.0x10140	7.0	10
CBA7.0x15140		15
CBA7.0x20140		20
CBA7.0x25140		25
CBA7.0x30140		30
CBA7.0x35140		35
CBA7.0x40140	40	

Stryderyk

CORONARY ANGIOPLASTY CATHETER WITH DELIVERY SYSTEM RAPID EXCHANGE (RX)



Technical data overview:

Balloon material	polyamide
Catheter coating	hydrophilic
Balloon nominal diameters	1.25 mm – 7.00 mm
Balloon nominal lengths	10 mm – 40 mm
Compatible guidewire	max. .014"
Compatible guiding catheter OD/ID	5F (.056")*
Inflation and deflation times	short
Tip entry profile	0.017" (0.43 mm)
Nominal pressure	10 atm**
Rated Burst Pressure (RBP)	16 atm***
Proximal shaft diameter	2F
Distal shaft diameter	2.7F
Catheter length	140 cm

*) for balloon diameter 5.0 mm, 6.0 mm and 7.0 mm use guiding catheter 6F

**) Nominal pressure 8 atm for balloon diameter 5.0 mm; 6.0 mm and 7.0 mm

***) Rated Burst Pressure (RBP) 12 atm for balloon diameter 5.0 mm; 6.0 mm and 7.0 mm

Caution:

The content of this document, including provided data, information, trademarks and logos are the sole property of BALTON. Any reproduction, whether in whole or in part, is strictly prohibited and will be considered a violation of BALTON's copyrights and other intellectual proprietary rights. ©2021 Balton Sp. z o.o. All rights reserved. This document and any images contained therein is intended solely for the use by healthcare professionals and BALTON's distributors.

This publication is for information purposes only and should not be construed as a commercial offer in the meaning of art. 66(1) of the Polish Civil Code.

All medical devices made commercially available by BALTON must be used in accordance with the enclosed instructions for use.

This publication presents only a small sample of our complete offering of medical devices. For further information please contact your Balton representative or email us at balton@balton.pl.



DISTRIBUIDOR EXCLUSIVO

Bioline Comercial Ltda.

R. José Getúlio, 579 14º andar cj 141/146

Liberdade | SP | 01509 001

+55 11 3274 7600

www.bioline.com.br