


RIVER[®] CTO

(for Chronic Total Occlusion)

CORONARY ANGIOPLASTY CATHETER
WITH DELIVERY SYSTEM RAPID EXCHANGE (RX)

World's first
smallest PTCA Balloon


Ø 0.75 mm



BALTON[®]



CE₂₄₀₉



**World's first smallest PTCA Balloon
for CTO procedure – 0.75 mm.**

FEATURES:

- Extraordinary flexibility
- Excellent pushability & trackability
- Unique tapered tip
- Two markers for precise positioning in the vessel

	Reference:	Nominal balloon diameter [mm]	Nominal balloon length [mm]	Radiopaque markers
MICRO BALLOONS	CBR0.75x10140CTO2	0.75	10	TWO MARKERS
	CBR0.75x15140CTO2		15	
	CBR0.75x20140CTO2		20	
	CBR0.75x25140CTO2		25	
	CBR0.75x30140CTO2		30	
	CBR0.75x35140CTO2		35	
	CBR0.75x40140CTO2		40	
	CBR0.80x10140CTO2	0.80	10	TWO MARKERS
	CBR0.80x15140CTO2		15	
	CBR0.80x20140CTO2		20	
	CBR0.80x25140CTO2		25	
	CBR0.80x30140CTO2		30	
CBR0.80x35140CTO2	35			
CBR0.80x40140CTO2	40			
CBR0.90x10140CTO2	0.90	10	TWO MARKERS	
CBR0.90x15140CTO2		15		
CBR0.90x20140CTO2		20		
CBR0.90x25140CTO2		25		
CBR0.90x30140CTO2		30		
CBR0.90x35140CTO2		35		
CBR0.90x40140CTO2		40		
STANDARD BALLOONS	CBR1.00x10140CTO2	1.00	10	TWO MARKERS
	CBR1.00x15140CTO2		15	
	CBR1.00x20140CTO2		20	
	CBR1.00x25140CTO2		25	
	CBR1.00x30140CTO2		30	
	CBR1.00x35140CTO2		35	
	CBR1.00x40140CTO2		40	
	CBR1.25x10140CTO2	1.25	10	TWO MARKERS
	CBR1.25x15140CTO2		15	
	CBR1.25x20140CTO2		20	
	CBR1.25x25140CTO2		25	
	CBR1.25x30140CTO2		30	
	CBR1.25x35140CTO2		35	
	CBR1.25x40140CTO2		40	

Minimum Guiding Catheter I.D. .056"

SINGLE MARKER CATHETERS ARE ALSO AVAILABLE

When ordering a single marker catheter please substitute **CTO2** with **CTO1** in the Ref. Number, e.g.: CBR0.75X30140**CTO1**.

Technical data overview:

Pushability and trackability	excellent
Inflation and deflation time	very short
Compatible guidewire	max. 0.14"
Compatible guiding catheter OD/ID	5F/0.056" (1.42 mm)
Entry profile	0.016"
Balloon nominal diameters	0.75-1.25 mm
Balloon nominal lengths	10-40 mm
Nominal pressure	8 atm
Rated Burst Pressure	18 atm
Catheter length	140 cm
Delivery system	Rapid Exchange (RX)

RIVER[®] CTO

(for Chronic Total Occlusion)

CORONARY ANGIOPLASTY CATHETER
WITH DELIVERY SYSTEM RAPID EXCHANGE (RX)

Other available balloon catheters:


Coronary angioplasty balloon catheters

RIVER[®]



Peripheral angioplasty balloon catheters

LOVIX[®]



Valvuloplasty balloon catheters

Valver[®]



Caution:

The content of this document, including provided data, information, trademarks and logos are the sole property of BALTON. Any reproduction, whether in whole or in part, is strictly prohibited and will be considered a violation of BALTON's copyrights and other intellectual proprietary rights. ©2021 Balton Sp. z o.o. All rights reserved. This document and any images contained therein is intended solely for the use by healthcare professionals and BALTON's distributors.

This publication is for information purposes only and should not be construed as a commercial offer in the meaning of art. 66(1) of the Polish Civil Code.

All medical devices made commercially available by BALTON must be used in accordance with the enclosed instructions for use.

This publication presents only a small sample of our complete offering of medical devices. For further information please contact your Balton representative or email us at balton@balton.pl.



Bioline Comercial Ltda.

R. José Getúlio, 579 14º andar cj 141/146

Liberdade | SP | 01509 001

+55 11 3274 7600

www.bioline.com.br