

# RIVER<sup>®</sup>

# RIVER<sup>®</sup>

*Coronary angioplasty balloon catheter*

**NEW, LOW PROFILE ANGIOPLASTY BALLOON CATHETER:**

- *Dedicated for narrow and tortuous vessels*
- *Extraordinary flexibility*
- *Excellent pushability and trackability*

Caution:

This publication could include technical or other inaccuracies or typographical errors. Changes are periodically added to the information herein; these changes will be incorporated in the new editions of the publication. Balton Ltd. may make improvements and/or changes in the product(s) described in this publication at any time.

01/2013




Bioline Comercial Ltda.  
R. José Getúlio, 579 14º andar cj 141/146  
Liberdade | SP | 01509 001  
+55 11 3274 7600  
[www.bioline.com.br](http://www.bioline.com.br)



# RIVER<sup>®</sup>

Coronary angioplasty balloon catheter



## Technical data :

Balloon material	nylon
Nominal pressure	8 atm
Catheter coating	hydrophilic
Tip entry profile	0.017" (0.43 mm)
Crossing profile	0.025" (0.64 mm)
Proximal shaft diameter	1.8F (0.60 mm)
Distal shaft diameter	2.5F (0.84 mm)
Pushability and trackability	excellent
Inflation and deflation time	very short
Compatible guidewire	max. 0.014"
Compatible guiding catheter OD/ID	5F / 0.056"
Catheter length	140 cm
Delivery system	Rapid Exchange

Reference:	Balloon diameter / length
CBR1.25x10140	1.25 / 10 mm
CBR1.25x15140	1.25 / 15 mm
CBR1.25x20140	1.25 / 20 mm
CBR1.25x25140	1.25 / 25 mm
CBR1.25x30140	1.25 / 30 mm
CBR1.25x35140	1.25 / 35 mm
CBR1.25x40140	1.25 / 40 mm
CBR1.50x10140	1.50 / 10 mm
CBR1.50x15140	1.50 / 15 mm
CBR1.50x20140	1.50 / 20 mm
CBR1.50x25140	1.50 / 25 mm
CBR1.50x30140	1.50 / 30 mm
CBR1.50x35140	1.50 / 35 mm
CBR1.50x40140	1.50 / 40 mm
CBR1.75x10140	1.75 / 10 mm
CBR1.75x15140	1.75 / 15 mm
CBR1.75x20140	1.75 / 20 mm
CBR1.75x25140	1.75 / 25 mm
CBR1.75x30140	1.75 / 30 mm
CBR1.75x35140	1.75 / 35 mm
CBR1.75x40140	1.75 / 40 mm
CBR2.00x10140	2.00 / 10 mm
CBR2.00x15140	2.00 / 15 mm
CBR2.00x20140	2.00 / 20 mm
CBR2.00x25140	2.00 / 25 mm
CBR2.00x30140	2.00 / 30 mm
CBR2.00x35140	2.00 / 35 mm
CBR2.00x40140	2.00 / 40 mm
CBR2.25x10140	2.25 / 10 mm
CBR2.25x15140	2.25 / 15 mm
CBR2.25x20140	2.25 / 20 mm
CBR2.25x25140	2.25 / 25 mm
CBR2.25x30140	2.25 / 30 mm
CBR2.25x35140	2.25 / 35 mm
CBR2.25x40140	2.25 / 40 mm
CBR2.50x10140	2.50 / 10 mm
CBR2.50x15140	2.50 / 15 mm
CBR2.50x20140	2.50 / 20 mm
CBR2.50x25140	2.50 / 25 mm
CBR2.50x30140	2.50 / 30 mm
CBR2.50x35140	2.50 / 35 mm
CBR2.50x40140	2.50 / 40 mm

Reference:	Balloon diameter / length
CBR2.75x10140	2.75 / 10 mm
CBR2.75x15140	2.75 / 15 mm
CBR2.75x20140	2.75 / 20 mm
CBR2.75x25140	2.75 / 25 mm
CBR2.75x30140	2.75 / 30 mm
CBR2.75x35140	2.75 / 35 mm
CBR2.75x40140	2.75 / 40 mm
CBR3.00x10140	3.00 / 10 mm
CBR3.00x15140	3.00 / 15 mm
CBR3.00x20140	3.00 / 20 mm
CBR3.00x25140	3.00 / 25 mm
CBR3.00x30140	3.00 / 30 mm
CBR3.00x35140	3.00 / 35 mm
CBR3.00x40140	3.00 / 40 mm
CBR3.25x10140	3.25 / 10 mm
CBR3.25x15140	3.25 / 15 mm
CBR3.25x20140	3.25 / 20 mm
CBR3.25x25140	3.25 / 25 mm
CBR3.25x30140	3.25 / 30 mm
CBR3.25x35140	3.25 / 35 mm
CBR3.25x40140	3.25 / 40 mm
CBR3.50x10140	3.50 / 10 mm
CBR3.50x15140	3.50 / 15 mm
CBR3.50x20140	3.50 / 20 mm
CBR3.50x25140	3.50 / 25 mm
CBR3.50x30140	3.50 / 30 mm
CBR3.50x35140	3.50 / 35 mm
CBR3.50x40140	3.50 / 40 mm
CBR3.75x10140	3.75 / 10 mm
CBR3.75x15140	3.75 / 15 mm
CBR3.75x20140	3.75 / 20 mm
CBR3.75x25140	3.75 / 25 mm
CBR3.75x30140	3.75 / 30 mm
CBR3.75x35140	3.75 / 35 mm
CBR3.75x40140	3.75 / 40 mm
CBR4.00x10140	4.00 / 10 mm
CBR4.00x15140	4.00 / 15 mm
CBR4.00x20140	4.00 / 20 mm
CBR4.00x25140	4.00 / 25 mm
CBR4.00x30140	4.00 / 30 mm
CBR4.00x35140	4.00 / 35 mm
CBR4.00x40140	4.00 / 40 mm